






Angle Rotors for CM-12 Plus

 Angle Rotor RM-1214 M (10 x 2ml) Max speed RPM 16000 Max RCF (g) 16600	
Image of Tube	
Tube Capacity (ml)	2
Type of Tube	Glass
Insert or Cushion No.	C 5R Included
No of tubes/Insert	1
No of tubes/rotor	10
Ø x L (mm)	12 x 36
Tubes supplied	Yes

 Angle Rotor RM-1215 M (24 x 1.5ml) Max speed RPM 11850 Max RCF (g) 12715	
Image of Tube	
Tube Capacity (ml)	1.5
Type of Tube	Micro
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	24
Ø x L (mm)	11 x 41
Tubes supplied	No

 Angle Rotor RM-1216 M (8 x 5ml) Max speed RPM 11850 Max RCF (g) 13500	
Image of Tube	
Tube Capacity (ml)	5
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	8
Ø x L (mm)	12 x 76
Tubes supplied	Yes

 Angle Rotor RM-1217 M (4 Strips of 8 PCR x 0.2ml) Max speed RPM 12000 Max RCF (g) 11910	
Image of Tube	
Tube Capacity (ml)	4 strips of 8 PCR
Type of Tube	PCR
No of tubes/rotor	4
Ø x L (mm)	6 x 78
Tubes supplied	No