

























Angle Rotors for CM-8 Plus

Angle Rotor R-83 M (16 x 15ml) Max. Speed RPM 6000 Max RCF (g) 4700 	Image of Tube						
	Tube Capacity (ml)	15	8	6	5	4	3
	Type of Tube	Plastic	Glass	Glass	Glass	Glass	Glass
	Insert or Cushion No.	C 15R Included	RE 15-8R (OLD - RA-154)	RE 15-6R (OLD - RA-156)	RE 15-5R (OLD - RA-153)	(OLD - RA-155)	RE 15-3R (OLD - RA-152)
	No of tubes/Insert	1	1	1	1		1
	No of tubes/rotor	16	16	16	16		16
	Ø x L (mm)	17 x 102	12 x 100	12 x 95	12 x 75	12 x 70	10 x 75
	Tubes supplied	Yes	Yes	No	Yes	No	Yes

Angle Rotor R-83 BM (8 x 15ml) Max. Speed RPM 6000 Max RCF (g) 4700 	Image of Tube						
	Tube Capacity (ml)	15	8	6	5	4	3
	Type of Tube	Plastic	Glass	Glass	Glass	Glass	Glass
	Insert or Cushion No.	C 15R Included	RE 15-8R (OLD - RA-154)	RE 15-6R (OLD - RA-156)	RE 15-5R (OLD - RA-153)	(OLD - RA-155)	RE 15-3R (OLD - RA-152)
	No of tubes/Insert	1	1	1	1		1
	No of tubes/rotor	8	8	8	8		8
	Ø x L (mm)	17 x 102	12 x 100	12 x 95	12 x 75	12 x 70	10 x 75
	Tubes supplied	Yes	Yes	No	Yes	No	Yes

Angle Rotor R-84 M (6 x 50ml) Max. Speed RPM 6000 Max RCF (g) 4700 	Image of Tube					
	Tube Capacity (ml)	50	30	15	15	10
	Type of Tube	Plastic	Glass	Glass	Conical	Blood Collection Tube
	Insert or Cushion No.	C 50R Included	(OLD - RA-502)	RE 50-15RR (OLD - RA-501)	RE 50-15FR (OLD - RA-501F)	RE 50-10R (OLD - RA-164)
	No of tubes/Insert	1		1	1	1
	No of tubes/rotor	6		6	6	6
	Ø x L (mm)	30 x 100	24 x 110	17 x 110	16.5 x 120	15.5 x 106
	Tubes supplied	Yes	Yes	Yes	No	No

Angle Rotor R-88 M (4 x 100ml) Max. Speed RPM 6000 Max RCF (g) 4550 	Image of Tube							
	Tube Capacity (ml)	100	50	50	15	15	15	10
	Type of Tube	Plastic	Conical	Glass	Glass	Glass	Conical	Blood Collection Tube
	Insert or Cushion No.	C 100R Included	RE 100-50FR (OLD - RA-161F)	RE 100-50RR (OLD - RA-161)	(OLD - RA-160)	RE 100-15R(3) (OLD - RA-166)	RE 100-15FR (OLD - RA-160F)	RE 100-10(3) (OLD - RA-168)
	No of tubes/Insert	1	1	1		3	1	3
	No of tubes/rotor	4	4	4		12	4	12
	Ø x L (mm)	45.5 x 104	29 x 116	30 x 110	17 x 110	17 x 110	16.5 x 120	15.5 x 106
	Tubes supplied	Yes	No	Yes	Yes	Yes	No	No

Angle Rotor R-89 M (24 x 15ml) Max. Speed RPM 6000 Max RCF (g) 5070 	Image of Tube						
	Tube Capacity (ml)	15	8	6	5	4	3
	Type of Tube	Plastic	Glass	Glass	Glass	Glass	Glass
	Insert or Cushion No.	C 15R Included	RE 15-8R (OLD - RA-154)	RE 15-6R (OLD - RA-156)	RE 15-5R (OLD - RA-153)	(OLD - RA-155)	RE 15-3R (OLD - RA-152)
	No of tubes/Insert	1	1	1	1		1
	No of tubes/rotor	24	24	24	24		24
	Ø x L (mm)	17 x 102	12 x 100	12 x 95	12 x 75	12 x 70	10 x 75
	Tubes supplied	Yes	Yes	No	Yes	No	Yes