


Angle Rotors for R-8C Plus & R-8M Plus



Angle Rotor A 4-100

Max. Speed RPM 6000

Max RCF (g) 4750

Image of Tube											
Tube Capacity (ml)	100	100	50	50	50	30	15	15	15	10	5/7
Type of Tube	Plastic	Nalgene® Oak Ridge	Conical	Glass	Nalgene® Oak Ridge	Nalgene® Oak Ridge	Glass	Glass	Conical	Blood Collection Tube	Blood Collection Tube (Long)
Insert or Cushion No.	C 100R Included	I 100NR	RE 100-50FR (OLD - RA-161F)	RE 100-50RR (OLD - RA-161)	RE 100-50NR	RE 100-30NR	(OLD - RA-160)	RE 100-15R(3) (OLD - RA-166)	RE 100-15FR (OLD - RA-160F)	RE 100-10(3) (OLD - RA-168)	RE 100-5/7L(4) (OLD - RA-167)
No of tubes/Insert	1	1	1	1	1	1	4	3	1	3	4
No of tubes/rotor	4	4	4	4	4	4	4	12	4	12	16
Ø x L (mm)	45.5 x 104	38 x 108.5	29 x 116	30 x 110	28 x 111	25 x 97	17 x 110	17 x 110	16.5 x 120	15.5 x 106	12.5 x 106
Tubes supplied	Yes	No	No	Yes	No	No	Yes	Yes	No	No	No

Angle Rotor AR 6-50

Max. Speed RPM 6000

Max RCF (g) 4870



Image of Tube							
Tube Capacity (ml)	50	50	50	30	15	15	10
Type of Tube	Plastic	Conical	Nalgene® Oak Ridge	Nalgene® Oak Ridge	Glass	Conical	Blood Collection Tube
Insert or Cushion No.	C 50R Included	C 50F	I 50N	RE 50-30NL	RE 50-15RR (OLD - RA-501)	RE 50-15FR (OLD-RA-501F)	RE 50-10R (OLD-RA-164)
No of tubes/Insert	1	1	1	1	1	1	1
No of tubes/rotor	6	6	6	6	6	6	6
Ø x L (mm)	30 x 100	29 x 116	28 x 111	25 x 97	17 x 110	16.5 x 120	15.5 x 106
Tubes supplied	Yes	No	No	No	Yes	No	No

Angle Rotor A 24-15

Max. Speed RPM 6000

Max RCF (g) 5070




Image of Tube							
Tube Capacity (ml)	15	15	8	5/7	5/7	5	3
Type of Tube	Plastic	Conical	Glass	Blood Collection Tube (Long)	Blood Collection Tube (Short)	Glass	Glass
Insert or Cushion No.	C 15R Included	C 15F	RE 15-8R (OLD - RA-154)	RE 15-5/7LR (OLD - RA-163)	RE 15-5/7SR (OLD - RA-162)	RE 15-5R (OLD - RA-153)	RE 15-3R (OLD - RA-152)
No of tubes/Insert	1	1	1	1	1	1	1
No of tubes/rotor	24	24	24	24	24	24	24
Ø x L (mm)	17 x 102	16.5 x 120	12 x 100	12.5 x 106	12.5 x 81	12 x 75	10 x 75
Tubes supplied	Yes	No	Yes	No	No	Yes	Yes

Angle Rotor A 12-15

Max. Speed RPM 6000

Max RCF (g) 5000




Image of Tube							
Tube Capacity (ml)	15	15	8	5/7	5/7	5	3
Type of Tube	Plastic	Conical	Glass	Blood Collection Tube (Long)	Blood Collection Tube (Short)	Glass	Glass
Insert or Cushion No.	C 15R Included	C 15F	RE 15-8R (OLD - RA-154)	RE 15-5/7LR (OLD - RA-163)	RE 15-5/7SR (OLD - RA-162)	RE 15-5R (OLD - RA-153)	RE 15-3R (OLD - RA-152)
No of tubes/Insert	1	1	1	1	1	1	1
No of tubes/rotor	12	12	12	12	12	12	12
Ø x L (mm)	17 x 102	16.5 x 120	12 x 100	12.5 x 106	12.5 x 81	12 x 75	10 x 75
Tubes supplied	Yes	No	Yes	No	No	Yes	Yes