





Swing Out Rotors for CM-8 Plus

Swing Out Rotor R-90 M (4 x 100ml)								
 Max. Speed RPM 4000 Max RCF (g) 2500				B-100				
Image of Tube								
Tube Capacity (ml)	100	50	50	30	15	15	10	5/7
Type of Tube	Glass	Glass	Conical	Glass	Glass	Conical	Blood Collection Tube	Blood Collection Tube (Long)
Insert or Cushion No.	C 100R Included	RE 100-50 (OLD - RA-157)	RE 100-50FR (OLD - RA-161F)	(OLD - RA-158)	RE 100-15R(3) (OLD - RA-166)	RE 100-15FR (OLD - RA-160F)	RE 100-10(3) (OLD - RA-168)	RE 100-5/7L(4) (OLD - RA-167)
No of tubes/Insert	1	1	1		3	1	3	4
No of tubes/rotor	4	4	4		12	4	12	16
Ø x L (mm)	45.5 x 104	30 x 110	29 x 116	24 x 110	17 x 110	16.5 x 120	15.5 x 106	12.5 x 106
Tubes supplied	Yes	Yes	No	Yes	Yes	No	No	No

Swing Out Rotor R-81 M (16 x 15ml)								
 Max. Speed RPM 4500 Max RCF (g) 3485								
Image of Tube								
Tube Capacity (ml)	15	8	5/7	5/7	6	5	4	3
Type of Tube	Glass	Glass	Blood Collection Tube (Long)	Blood Collection Tube (Short)	Glass	Glass	Glass	Wintrob
Insert or Cushion No.	C 15R Included	RE 15-8R (OLD - RA-154)	RE 15-5/7LR (OLD - RA-163)	RE 15-5/7SR (OLD - RA-162)	RE 15-6R (OLD - RA-156)	RE 15-5R (OLD - RA-153)	(OLD - RA-155)	RE 15-3R (OLD - RA-152)
No of tubes/Insert	1	1	1	1	1	1		1
No of tubes/rotor	16	16	16	16	16	16	16	16
Ø x L (mm)	17 x 110	12 x 100	12.5 x 106	12.5 x 81	12 x 95	12 x 75	10 x 70	10 x 75
Tubes supplied	Yes	Yes	No	No	No	Yes	No	Yes

Swing Out Rotor R-81 BM (8 x 15ml)								
 Max. Speed RPM 4500 Max RCF (g) 3485								
Image of Tube								
Tube Capacity (ml)	15	8	5/7	5/7	6	5	4	3
Type of Tube	Glass	Glass	Blood Collection Tube (Long)	Blood Collection Tube (Short)	Glass	Glass	Glass	Wintrob
Insert or Cushion No.	C 15R Included	RE 15-8R (OLD - RA-154)	RE 15-5/7LR (OLD - RA-163)	RE 15-5/7SR (OLD - RA-162)	RE 15-6R (OLD - RA-156)	RE 15-5R (OLD - RA-153)	(OLD - RA-155)	RE 15-3R (OLD - RA-152)
No of tubes/Insert	1	1	1	1	1	1		1
No of tubes/rotor	8	8	8	8	8	8	8	8
Ø x L (mm)	17 x 110	12 x 100	12.5 x 106	12.5 x 81	12 x 95	12 x 75	10 x 70	10 x 75
Tubes supplied	Yes	Yes	No	No	No	Yes	No	Yes



Swing Out Rotor R-82 AM (4 x 50ml)

Max. Speed RPM 4500

Max RCF (g) 3485

Image of Tube					
Tube Capacity (ml)	50	30	15	15	10
Type of Tube	Glass	Glass	Glass	Conical	Blood Collection Tube
Insert or Cushion No.	C 15R Included	(OLD - RA-502)	RE 50-15RR (OLD - RA-501)	RE 50-15FR (OLD - RA-501F)	RE 50-10R (OLD - RA-164)
No of tubes/Insert	1	1	1	1	1
No of tubes/rotor	4	4	4	4	4
Ø x L (mm)	30 x 110	24 x 110	17 x 110	16.5 x 120	15.5 x 106
Tubes supplied	Yes	Yes	Yes	No	No

Discontinued