



Rotors for RM-03 Pro

 Micro Tube Rotor Max speed RPM 14500 Max RCF (g) 14104				 4 Strips PCR Rotor Max. Speed RPM 14500 Max RCF (g) 12576	
Image of Tube				Image of Tube	
Tube Capacity (ml)	2	0.5	0.2	Tube Capacity (ml)	4 Strips of 8 x 0.2
Type of Tube	Micro	Micro	Micro	Type of Tube	PCR
Insert Supplied	No	Yes	Yes	No of tubes/rotor	4
No of tubes/insert	-	1	1	Ø x L (mm)	6 x 78
No of tubes/rotor	12	12	12	Tubes supplied	No
Ø x L (mm)	11 x 41	7.5 x 31	6 x 22		
Tubes supplied	No	No	No		

(Rotor supplied with all suitable inserts)