


Rotor for RM-02 Pro



Micro Hybrid Rotor

Max speed RPM 7000

Max RCF (g) 3280

Image of Tube				
Tube Capacity (ml)	2 / 1.5	0.5	0.2	4 Strips of 8 x 0.2
Type of Tube	Micro	Micro	Micro	PCR
Insert Supplied	No	Yes	Yes	No
No of tubes/Insert	-	1	1	-
No of tubes/rotor	12	12	12	4
Ø x L (mm)	11 x 41	7.5 x 31	6 x 22	6 x 78
Tubes supplied	No	No	No	No

(Rotor supplied with all suitable inserts)