















Rotors Listing for R-23 & PR-23


 R-230		 R-231			Swing Out Rotor R-230 + R-231 (4 x 200ml)			
					R-23		PR-23	
					Max. Speed RPM 4800 Max RCF (g) 4000		Max. Speed RPM 5100 Max RCF (g) 4500	
Image of Tube								
Tube Capacity (ml)	200	100	50	30	15	10	5/7	5/7
Type of Tube	Glass	Glass	Glass	Glass	Glass	Blood Collection Tube	Blood Collection Tube (Long)	Blood Collection Tube (Short)
Insert or Cushion No.	C 200R Included	(OLD - RA-2311)	(OLD - RA-2312)	(OLD - RA-2313)	RE 200-15R(7) (OLD - RA-2314)	(OLD - RA-2316)	(OLD - RA-2315)	(OLD - RA-2317)
Image of Insert (Set of 4, supplied with suitable bottom cushions)	No							
No of tubes/Insert	1	1	1	3	7	5	6	6
No of tubes/rotor	4	4	4	12	28	20	24	24
Ø x L (mm)	55 x 115	45.5 x 104	30 x 110	24 x 110	17 x 110	15.5 x 106	12.5 x 106	12.5 x 81
Tubes supplied	Yes	Yes	Yes	Yes	Yes	No	No	No


Swing Out Rotor R-230B + R-233B (108 X 5 (RIA))	
R-23	PR-23
Max. Speed RPM 5000 Max RCF (g) 3830	Max. Speed RPM 5100 Max RCF (g) 3980
	
Image of Tube	
Tube Capacity (ml)	108 X 5 (RIA)
Type of Tube	RIA Plastic Tubes
Insert or Cushion No.	No
No of tubes/Carrier	27
No of tubes/rotor	108
Ø x L (mm)	12.4 x 76
Tubes supplied	Yes


		Swing Out Rotor R-234 (4 x 250ml)			
		PR-23			
		Max. Speed RPM 5500 Max RCF (g) 5070			
Image of Tube					
Tube Capacity (ml)	250	50	25	15	5
Type of Tube	Plastic	Glass	Glass	Glass	Glass
Insert or Cushion No.	C 250R Included	(OLD - RA-2341)	(OLD - RA-2342)	RE 250-15R(5) (OLD - RA-2343)	RE 250-5R(13) (OLD - RA-2344)
Image of Insert (Set of 4, supplied with suitable bottom cushions)	No				
No of tubes/Insert	1	1	1	5	13
No of tubes/rotor	4	4	2	20	52
Ø x L (mm)	67.5 x 100	30 x 110	28 x 62.5	17 x 110	14 x 100
Tubes supplied	Yes	Yes	Yes	Yes	Yes



Swing Out Rotor R-235 (8 x 96 WELL)	
	PR-23
	Max. Speed RPM 3200 Max RCF (g) 1470
Image of Tube	
Tube Capacity (ml)	8 x 96 wells
Type of Tube	Microtiter
No of tubes/rotor	8
W x D x L (mm)	127 x 84 x 18
Tubes supplied	No


Swing Out Rotor R-230A (2 x 96 WELL)	
	PR-23
	Max. Speed RPM 3200 Max RCF (g) 1470
Image of Tube	
Tube Capacity (ml)	2 x 96 wells
Type of Tube	Microtiter
No of tubes/rotor	2
W x D x L (mm)	127 x 84 x 18
Tubes supplied	No



Angle Rotor R-236 (6 x 200ml)	
R-23	PR-23
Max. Speed RPM 6900 Max RCF (g) 7350	Max. Speed RPM 7000 Max RCF (g) 7560
	
Image of Tube	
Tube Capacity (ml)	200
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	6
Ø x L (mm)	61 x 125
Tubes supplied	Yes

Angle Rotor R-237 (8 x 100ml)	
R-23	PR-23
Max. Speed RPM 6900 Max RCF (g) 7200	Max. Speed RPM 7000 Max RCF (g) 7400
	
Image of Tube	
Tube Capacity (ml)	100
Type of Tube	Plastic
Insert or Cushion No.	C 100R Included
No of tubes/Insert	1
No of tubes/rotor	8
Ø x L (mm)	45.5 x 104
Tubes supplied	Yes

Angle Rotor R-238 (12 x 50ml)	
R-23	PR-23
Max. Speed RPM 6900 Max RCF (g) 6970	Max. Speed RPM 7000 Max RCF (g) 7450
	
Image of Tube	
Tube Capacity (ml)	50
Type of Tube	Plastic
Insert or Cushion No.	C 50R Included
No of tubes/Insert	1
No of tubes/rotor	12
Ø x L (mm)	30 x 100
Tubes supplied	Yes

Angle Rotor R-239 (36 x 15ml)	
R-23	PR-23
Max. Speed RPM 6900 Max RCF (g) 7240	Max. Speed RPM 7000 Max RCF (g) 7400
	
Image of Tube	
Tube Capacity (ml)	15
Type of Tube	Plastic
Insert or Cushion No.	C 15R Included
No of tubes/Insert	1
No of tubes/rotor	36
Ø x L (mm)	17 x 102
Tubes supplied	Yes

Angle Rotor R-238F (12 x 50ml)	
R-23	PR-23
Max. Speed RPM 6900 Max RCF (g) 6970	Max. Speed RPM 7000 Max RCF (g) 7450
	
Image of Tube	
Tube Capacity (ml)	50
Type of Tube	Conical
Insert or Cushion No.	C 50F
No of tubes/Insert	1
No of tubes/rotor	12
Ø x L (mm)	29 x 116
Tubes supplied	No

Angle Rotor R-239F (24 x 15ml)	
R-23	PR-23
Max. Speed RPM 6900 Max RCF (g) 7240	Max. Speed RPM 7000 Max RCF (g) 7400
	
Image of Tube	
Tube Capacity (ml)	15
Type of Tube	Conical
Insert or Cushion No.	C 15F
No of tubes/Insert	1
No of tubes/rotor	24
Ø x L (mm)	16.5 x 120
Tubes supplied	No