






Rotors Listing for R-24 & PR-24


Swing Out Rotor R-246A (4 x 25ml)	
R-24	PR-24
Max. Speed RPM 9300 Max RCF (g) 8700	Max. Speed RPM 11000 Max RCF (g) 12180
	
Image of Tube	
Tube Capacity (ml)	25
Type of Tube	Glass
Insert or Cushion No.	C 25R Included
No of tubes/Insert	1
No of tubes/rotor	4
Ø x L (mm)	28 x 62.5
Tubes supplied	Yes



Angle Rotor R-241 (4 x 100ml)	
R-24	PR-24
Max. Speed RPM 11000 Max RCF (g) 14050	Max. Speed RPM 14500 Max RCF (g) 23980
	
Image of Tube	
Tube Capacity (ml)	100
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	4
Ø x L (mm)	45.5 x 104
Tubes supplied	Yes


Angle Rotor R-242 (6 x 50ml)	
R-24	PR-24
Max. Speed RPM 12100 Max RCF (g) 14240	Max. Speed RPM 16000 Max RCF (g) 25180
	
Image of Tube	
Tube Capacity (ml)	50
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	6
Ø x L (mm)	30 x 100
Tubes supplied	Yes


Angle Rotor R-243 (8 x 25ml)	
R-24	PR-24
Max. Speed RPM 12500 Max RCF (g) 14500	Max. Speed RPM 14500 Max RCF (g) 20220
	
Image of Tube	
Tube Capacity (ml)	25
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	8
Ø x L (mm)	25 x 95
Tubes supplied	Yes



Angle Rotor R-244 (12 x 15ml)	
R-24	PR-24
Max. Speed RPM 12300 Max RCF (g) 14710	Max. Speed RPM 16000 Max RCF (g) 24900
	
Image of Tube	
Tube Capacity (ml)	15
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	12
Ø x L (mm)	17 x 102
Tubes supplied	Yes


Angle Rotor R-243A (6 x 25ml)	
R-24	PR-24
Max. Speed RPM 14100 Max RCF (g) 17120	Max. Speed RPM 18000 Max RCF (g) 28980
	
Image of Tube	
Tube Capacity (ml)	25
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	6
Ø x L (mm)	25 x 95
Tubes supplied	Yes


Angle Rotor R-244A (8 x 15ml)	
R-24	PR-24
Max. Speed RPM 15400 Max RCF (g) 20150	Max. Speed RPM 18000 Max RCF (g) 27530
	
Image of Tube	
Tube Capacity (ml)	15
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	8
Ø x L (mm)	17 x 102
Tubes supplied	Yes


Angle Rotor R-247 (8 x 50ml)	
R-24	PR-24
Max. Speed RPM 10800 Max RCF (g) 13430	Max. Speed RPM 14500 Max RCF (g) 24680
	
Image of Tube	
Tube Capacity (ml)	50
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	8
Ø x L (mm)	30 x 100
Tubes supplied	Yes

Angle Rotor R-248 (24 x 2ml)	
R-24	PR-24
Max. Speed RPM 17300 Max RCF (g) 27440	Max. Speed RPM 20000 Max RCF (g) 37570
	
Image of Tube	
Tube Capacity (ml)	2
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	24
Ø x L (mm)	11 x 41
Tubes supplied	Yes

Angle Rotor R-249 (6 Strips of 8 PCR x 0.2ml)	
R-24	PR-24
Max. Speed RPM 11000 Max RCF (g) 12180	Max. Speed RPM 11000 Max RCF (g) 12180
	
Image of Tube	
Tube Capacity (ml)	6 strips of 8 PCR
Type of Tube	PCR
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	6
Ø x L (mm)	6 x 78
Tubes supplied	Yes

Angle Rotor R-251 (28 x 1.5ml)	
R-24	PR-24
Max. Speed RPM 17300 Max RCF (g) 30300	Max. Speed RPM 20000 Max RCF (g) 40500
	
Image of Tube	
Tube Capacity (ml)	1.5
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	28
Ø x L (mm)	11 x 41
Tubes supplied	Yes

Angle Rotor R-244F (12 x 15ml)	
R-24	PR-24
Max. Speed RPM 11000 Max RCF (g) 15200	Max. Speed RPM 11000 Max RCF (g) 15200
	
Image of Tube	
Tube Capacity (ml)	15
Type of Tube	Conical
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	12
Ø x L (mm)	16.5 x 120
Tubes supplied	No

Angle Rotor R-247F (8 x 50ml)	
R-24	PR-24
Max. Speed RPM 11000 Max RCF (g) 15200	Max. Speed RPM 11000 Max RCF (g) 15200
	
Image of Tube	
Tube Capacity (ml)	50
Type of Tube	Conical
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	8
Ø x L (mm)	29 x 116
Tubes supplied	No