





Swing Out Rotors for R-8C Plus & R-8M Plus

Image of Tube	 Swing Out Rotor S 4-100 Max. Speed RPM 4000 Max RCF (g) 2500				Metal Buckets suitable for Swing Out Rotor S 4-100						
					B 100			B 9-5			
					 						
Tube Capacity (ml)	100	100	50	50	50	30	15	15	10	5/7	5/7
Type of Tube	Glass	Nalgene® Oak Ridge	Glass	Conical	Nalgene® Oak Ridge	Nalgene® Oak Ridge	Glass	Conical	Blood Collection Tube	Blood Collection Tube (Long)	Blood Collection Tube (Short)
Bucket No.	B 100	B 100	B 100	B 100	B 100	B 100	B 100	B 100	B 100	B 100	B 9-5
Insert or Cushion No.	C 100R Included	I 100N	RE 100-50 (OLD - RA-157)	RE 100-50FR (OLD - RA-161F)	RE 100-50N	RE 100-30N	RE 100-15R(3) (OLD - RA-166)	RE 100-15FR (OLD - RA-160F)	RE 100-10(3) (OLD - RA-168)	RE 100-5/7L(4) (OLD - RA-167)	-
No of tubes/Insert	1	1	1	1	1	1	3	1	3	4	9
No of tubes/rotor	4	4	4	4	4	4	12	4	12	16	36
Ø x L (mm)	45.5 x 104	38 x 108.5	30 x 110	29 x 116	28 x 111	25 x 97	17 x 110	16.5 x 120	15.5 x 106	12.5 x 106	12.5 x 81
Tubes supplied	Yes	No	Yes	No	No	No	Yes	No	No	No	No

Swing Out Rotor S4-100 with max. speed of 4500rpm, Max rcf (g) 3170 also available on Request

Swing Out Rotor S 8-15 Max. Speed RPM 4500 Max RCF (g) 3440 	Image of Tube							
	Tube Capacity (ml)	15	15	8	5/7	5/7	5	3
	Type of Tube	Glass	Conical	Glass	Blood Collection Tube(Long)	Blood Collection Tube (Short)	Glass	Glass
	Insert or Cushion No.	C 15R Included	C 15F	RE 15-8R (OLD - RA-154)	RE 15-5/7LR (OLD - RA-163)	RE 15-5/7SR (OLD - RA-162)	RE 15-5R (OLD - RA-153)	RE 15-3R (OLD - RA-152)
	No of tubes/Insert	1	1	1	1	1	1	1
	No of tubes/rotor	8	8	8	8	8	8	8
	Ø x L (mm)	17 x 110	16.5 x 120	12 x 100	12.5 x 106	12.5 x 81	12 x 75	10 x 75
	Tubes supplied	Yes	No	Yes	No	No	Yes	Yes

Swing Out Rotor S 6-50 Max. Speed RPM 4500 Max RCF (g) 3530 	Image of Tube							
	Tube Capacity (ml)	50	50	50	30	15	15	10
	Type of Tube	Glass	Conical	Nalgene® Oak Ridge	Nalgene® Oak Ridge	Conical	Glass	Blood Collection Tube
	Insert or Cushion No.	C 50R Included	C 50F	I 50N	RE 50-30NS	RE 50-15FR (OLD - RA-501F)	RE 50-15RR (OLD - RA-501)	RE 50-10R (OLD - RA-164)
	No of tubes/Insert	1	1	1	1	1	1	1
	No of tubes/rotor	6	6	6	6	6	6	6
	Ø x L (mm)	30 x 110	29 x 116	28 x 111	25 x 97	16.5 x 120	17 x 110	15.5 x 106
	Tubes supplied	Yes	No	No	No	No	Yes	No

Swing Out Rotor S 24-15 Max. Speed RPM 4500 Max RCF (g) 3330 	Image of Tube						
	Tube Capacity (ml)	15	8	5/7	5/7	5	3
	Type of Tube	Glass	Glass	Blood Collection Tube(Long)	Blood Collection Tube (Short)	Glass	Glass
	Insert or Cushion No.	C 15R Included	RE 15-8R (OLD - RA-154)	RE 15-5/7LR (OLD - RA-163)	RE 15-5/7SR (OLD - RA-162)	RE 15-5R (OLD - RA-153)	RE 15-3R (OLD - RA-152)
	No of tubes/Insert	1	1	1	1	1	1
	No of tubes/rotor	24	24	24	24	24	24
	Ø x L (mm)	17 x 110	12 x 100	12.5 x 106	12.5 x 81	12 x 75	10 x 75
	Tubes supplied	Yes	Yes	No	No	Yes	Yes

Specifications subject to change without notice

Swing Out Rotor S 16-15

Max. Speed RPM 4500

Max RCF (g) 3400




Image of Tube							
Tube Capacity (ml)	15	15	8	5/7	5/7	5	3
Type of Tube	Glass	Conical	Glass	Blood Collection Tube (Long)	Blood Collection Tube (Short)	Glass	Glass
Insert or Cushion No.	C 15R Included	C 15F	RE 15-8R (OLD – RA-154)	RE 15-5/7LR (OLD – RA-163)	RE 15-5/7SR (OLD – RA-162)	RE 15-5R (OLD – RA-153)	RE 15-3R (OLD – RA-152)
No of tubes/Insert	1	1	1	1	1	1	1
No of tubes/rotor	16	16	16	16	16	16	16
Ø x L (mm)	17 x 110	16.5 x 120	12 x 100	12.5 x 106	12.5 x 81	12 x 75	10 x 75
Tubes supplied	Yes	No	Yes	No	No	Yes	Yes