

Rotor Listing				Reduction Adaptors (With Glass Tubes)					
Models	Capacity	Max. Speed	max. RCF	RE 50-15R (OLD RA-501)	RA-50	RA-151	RE 50-15R (OLD RA-152)	RE 50-15R (OLD RA-153)	RE 50-15R (OLD RA-154)
		RPM	'g'	(1X15ml)	(1X30ml)	(1Xwintrob)	(1X3ml)	(1X5ml)	(1X8ml)
Swing out Heads (with glass tubes)									
R-41	8 X 15ml	3500	1850	X	X	X	✓	✓	✓
R-42	4 X 50ml	3500	1850	✓	✓	X	X	X	X
Angle Heads (with polypropylene tubes)									
R-43	8 X 15ml	4400	2200	X	X	✓	✓	✓	✓
R-44	4 X 50ml	4200	2350	✓	✓	X	X	X	X

Specifications subject to change without notice