






Angle Rotors


Angle Rotor R-236 M (6 x 200ml)	
C-23 Plus	C-30 Plus CPR-30 Plus
Max. Speed RPM 7000 Max RCF (g) 7400	Max. Speed RPM 7000 Max RCF (g) 7400
	
Image of Tube	
Tube Capacity (ml)	200
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	6
Ø x L (mm)	61 x 125
Tubes supplied	Yes


Angle Rotor R-237 M (8 x 100ml)	
C-23 Plus	C-30 Plus CPR-30 Plus
Max. Speed RPM 7000 Max RCF (g) 7400	Max. Speed RPM 7000 Max RCF (g) 7400
	
Image of Tube	
Tube Capacity (ml)	100
Type of Tube	Plastic
Insert or Cushion No.	C 100R Included
No of tubes/Insert	1
No of tubes/rotor	8
Ø x L (mm)	45.5 x 104
Tubes supplied	Yes

Angle Rotor R-238 M (12 x 50ml)	
C-23 Plus	C-30 Plus CPR-30 Plus
Max. Speed RPM 7000 Max RCF (g) 7400	Max. Speed RPM 7000 Max RCF (g) 7400
	
Image of Tube	
Tube Capacity (ml)	50
Type of Tube	Plastic
Insert or Cushion No.	C 50R Included
No of tubes/Insert	1
No of tubes/rotor	12
Ø x L (mm)	30 x 100
Tubes supplied	Yes


Angle Rotor R-239 M (36 x 15ml)	
C-23 Plus	C-30 Plus CPR-30 Plus
Max. Speed RPM 7000 Max RCF (g) 7400	Max. Speed RPM 7000 Max RCF (g) 7400
	
Image of Tube	
Tube Capacity (ml)	15
Type of Tube	Plastic
Insert or Cushion No.	C 15R Included
No of tubes/Insert	1
No of tubes/rotor	36
Ø x L (mm)	17 x 102
Tubes supplied	Yes

Angle Rotor R-241 M (4 x 100ml)	
C-24 Plus CPR-24 Plus	C-30 Plus CPR-30 Plus
Max. Speed RPM 14500 Max RCF (g) 23980	
	
Image of Tube	
Tube Capacity (ml)	100
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	4
Ø x L (mm)	45.5 x 104
Tubes supplied	Yes

Angle Rotor R-243 M (8 x 25ml)	
C-24 Plus CPR-24 Plus	C-30 Plus CPR-30 Plus
Max. Speed RPM 18000 Max RCF (g) 31150	
	
Image of Tube	
Tube Capacity (ml)	25
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	8
Ø x L (mm)	25 x 95
Tubes supplied	Yes


Angle Rotor R-244 M (12 x 15ml)	
C-24 Plus CPR-24 Plus	C-30 Plus CPR-30 Plus
Max. Speed RPM 18000 Max RCF (g) 31150	
	
Image of Tube	
Tube Capacity (ml)	15
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	12
Ø x L (mm)	17 x 102
Tubes supplied	Yes


Angle Rotor R-247 M (8 x 50ml)	
C-24 Plus CPR-24 Plus	C-30 Plus CPR-30 Plus
Max. Speed RPM 16000 Max RCF (g) 29480	
	
Image of Tube	
Tube Capacity (ml)	50
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	8
Ø x L (mm)	30 x 100
Tubes supplied	Yes


Angle Rotor R-248 M (24 x 2ml)	
C-24 Plus CPR-24 Plus	C-30 Plus CPR-30 Plus
Max. Speed RPM 20000 Max RCF (g) 37570	
	
Image of Tube	
Tube Capacity (ml)	2
Type of Tube	Plastic
Insert or Cushion No.	No
No of tubes/Insert	-
No of tubes/rotor	24
Ø x L (mm)	11 x 41
Tubes supplied	Yes


Specifications subject to change without notice.


Specifications subject to change without notice.

Angle Rotor R-249 M (6 Strips of 8 PCR x 0.2ml)		
C-24 Plus CPR-24 Plus	CPR-24i Plus	C-30 Plus CPR-30 Plus
 Max. Speed RPM 11000 Max RCF (g) 12180		
Image of Tube		
Tube Capacity (ml)	6 strips of 8 PCR	
Type of Tube	PCR	
Insert or Cushion No.	No	
No of tubes/Insert	-	
No of tubes/rotor	6	
Ø x L (mm)	6 x 78	
Tubes supplied	Yes	


Angle Rotor R-250 M (48 x 1.5ml)		
C-24 Plus CPR-24 Plus	CPR-24i Plus	C-30 Plus CPR-30 Plus
 Max. Speed RPM 16000 Max RCF (g) 27480		
Image of Tube		
Tube Capacity (ml)	1.5	
Type of Tube	Plastic	
Insert or Cushion No.	No	
No of tubes/Insert	-	
No of tubes/rotor	48	
Ø x L (mm)	11 x 41	
Tubes supplied	Yes	

Angle Rotor R-251 M (28 x 1.5ml)		
C-24 Plus CPR-24 Plus	CPR-24i Plus	C-30 Plus CPR-30 Plus
 Max. Speed RPM 20000 Max RCF (g) 37570		
Image of Tube		
Tube Capacity (ml)	1.5	
Type of Tube	Plastic	
Insert or Cushion No.	No	
No of tubes/Insert	-	
No of tubes/rotor	28	
Ø x L (mm)	11 x 41	
Tubes supplied	Yes	

Angle Rotor R-244 MF (12 x 15ml)		
C-24 Plus CPR-24 Plus	CPR-24i Plus	C-30 Plus CPR-30 Plus
 Max. Speed RPM 14500 Max RCF (g) 28680		
Image of Tube		
Tube Capacity (ml)	15	
Type of Tube	Conical	
Insert or Cushion No.	No	
No of tubes/Insert	-	
No of tubes/rotor	12	
Ø x L (mm)	16.5 x 120	
Tubes supplied	No	

Angle Rotor R-247 MF (8 x 50ml)		
C-24 Plus CPR-24 Plus	CPR-24i Plus	C-30 Plus CPR-30 Plus
 Max. Speed RPM 14500 Max RCF (g) 28680		
Image of Tube		
Tube Capacity (ml)	50	
Type of Tube	Conical	
Insert or Cushion No.	No	
No of tubes/Insert	-	
No of tubes/rotor	8	
Ø x L (mm)	29 x 116	
Tubes supplied	No	

Angle Rotor R-241 M (4 x 100ml)		
C-24 Plus CPR-24 Plus	CPR-24i Plus	C-30 Plus CPR-30 Plus
 Max. Speed RPM 14500 Max RCF (g) 23980		
Image of Tube		
Tube Capacity (ml)	100	
Type of Tube	Plastic	
Insert or Cushion No.	No	
No of tubes/Insert	-	
No of tubes/rotor	4	
Ø x L (mm)	45.5 x 104	
Tubes supplied	Yes	

Angle Rotor R-243 M (8 x 25ml)		
C-24 Plus CPR-24 Plus	CPR-24i Plus	C-30 Plus CPR-30 Plus
 Max. Speed RPM 18000 Max RCF (g) 31150		
Image of Tube		
Tube Capacity (ml)	25	
Type of Tube	Plastic	
Insert or Cushion No.	No	
No of tubes/Insert	-	
No of tubes/rotor	8	
Ø x L (mm)	30 x 95	
Tubes supplied	Yes	