

Models	C-23 Plus	C-24 Plus CPR-24 Plus	CPR-24i Plus	C-30 Plus CPR-30 Plus
 R-230 M + R-231 M 4 x 200 ml	Max 5100 rpm Max 4500 g	×	×	Max 4500 rpm Max 3500 g
 R-234 M 4 x 250 ml	Max 5500 rpm Max 5070 g	×	Max 4500 rpm Max 3500 g	Max 4500 rpm Max 3500 g
 R-235 M 8 x 96 wells	Max 3200 rpm Max 1370 g	×	Max 3200 rpm Max 1370 g	Max 3200 rpm Max 1370 g
 R-230 BM + R 233 BM 108 x 5 ml (RIA)	Max 5100 rpm Max 3980 g	×	×	Max 4500 rpm Max 3500 g
 R-240 M 16 x 50 ml	×	×	×	Max 3200 rpm Max 1775 g
 R-240 AM 28 x 15 ml	×	×	×	Max 3200 rpm Max 2110 g
 R-246 AM 4 x 25 ml	×	Max 11000 rpm Max 12180 g	Max 11000 rpm Max 12180 g	Max 11000 rpm Max 12180 g
 R-236 M 6 x 200 ml	Max 7000 rpm Max 7400 g	×	×	Max 7000 rpm Max 7400 g
 R-237 M 8 x 100 ml	Max 7000 rpm Max 7400 g	×	×	Max 7000 rpm Max 7400 g
 R-238 M 12 x 50 ml	Max 7000 rpm Max 7400 g	×	×	Max 7000 rpm Max 7400 g
 R-239 M 36 x 15 ml	Max 7000 rpm Max 7400 g	×	×	Max 7000 rpm Max 7400 g
 R-241 M 4 x 100 ml	×	Max 14500 rpm Max 23980 g	Max 14500 rpm Max 23980 g	Max 14500 rpm Max 23980 g
 R-243 M 8 x 25 ml	×	Max 18000 rpm Max 31150 g	Max 18000 rpm Max 31150 g	Max 18000 rpm Max 31150 g
 R-244 M 12 x 15 ml	×	Max 18000 rpm Max 31150 g	Max 18000 rpm Max 31150 g	Max 18000 rpm Max 31150 g
 R-247 M 8 x 50 ml	×	Max 16000 rpm Max 29480 g	Max 16000 rpm Max 29480 g	Max 16000 rpm Max 29480 g
 R-248 M 24 x 2 ml	×	Max 20000 rpm Max 37570 g	Max 20000 rpm Max 37570 g	Max 20000 rpm Max 37570 g
 R-249 M 6 Strips of 8 PCR	×	Max 11000 rpm Max 12180 g	Max 11000 rpm Max 12180 g	Max 11000 rpm Max 12180 g
 R-250 M 48 x 1.5 ml	×	Max 16000 rpm Max 30020 g	Max 16000 rpm Max 30020 g	Max 16000 rpm Max 30020 g
 R-251 M 28 x 1.5 ml	×	Max 20000 rpm Max 37570 g	Max 20000 rpm Max 37570 g	Max 20000 rpm Max 37570 g
 R-238 MF 12 x 50 ml	Max 7000 rpm Max 7400 g	×	×	Max 7000 rpm Max 7400 g
 R-239 MF 24 x 15 ml	Max 7000 rpm Max 7400 g	×	×	Max 7000 rpm Max 7400 g
 R-244 MF 12 x 15 ml	×	Max 14500 rpm Max 28680 g	Max 14500 rpm Max 28680 g	Max 14500 rpm Max 28680 g
 R-247 MF 8 x 50 ml	×	Max 14500 rpm Max 28680 g	Max 14500 rpm Max 28680 g	Max 14500 rpm Max 28680 g

Specifications subject to change without notice