












Swing Out Rotors



 R-230 M		 R-231 M		Swing Out Rotor R-230 M + R-231 M (4 x 200ml)				
				C-23 Plus		C-30 Plus CPR-30 Plus		
				Max. Speed RPM 5100 Max RCF (g) 4500		Max. Speed RPM 4500 Max RCF (g) 3500		
Image of Tube								
Tube Capacity (ml)	200	100	50	30	15	10	5/7	5/7
Type of Tube	Glass	Glass	Glass	Glass	Glass	Blood Collection Tube	Blood Collection Tube (Long)	Blood Collection Tube (Short)
Insert or Cushion No.	C 200R Included	(OLD - RA-2311)	(OLD - RA-2312)	(OLD - RA-2313)	RE 200-15R(7) (OLD - RA-2314)	(OLD - RA-2316)	(OLD - RA-2315)	(OLD - RA-2317)
Image of Insert (Set of 4, supplied with suitable bottom cushions)	No							
No of tubes/Insert	1	1	1	3	7	5	6	6
No of tubes/rotor	4	4	4	12	28	20	24	24
Ø x L (mm)	55 x 115	45.5 x 104	30 x 110	24 x 110	17 x 110	15.5 x 106	12.5 x 106	12.5 x 81
Tubes supplied	Yes	Yes	Yes	Yes	Yes	No	No	No


Swing Out Rotor R-234 M (4 x 250ml)					
C-23 Plus		CPR-24i Plus		C-30 Plus CPR-30 Plus	
Max. Speed RPM 5500 Max RCF (g) 5070				Max. Speed RPM 4500 Max RCF (g) 3500	
Image of Tube					
Tube Capacity (ml)	250	50	25	15	5
Type of Tube	Plastic	Glass	Glass	Glass	Glass
Insert or Cushion No.	C 250R Included	(OLD - RA-2341)	(OLD - RA-2342)	RE 250-15R(5) (OLD - RA-2343)	RE 250-5R(13) (OLD - RA-2344)
Image of Insert (Set of 4, supplied with suitable bottom cushions)	No				
No of tubes/Insert	1	1	1	5	13
No of tubes/rotor	4	4	4	20	52
Ø x L (mm)	67.5 x 100	30 x 110	25 x 62.5	17 x 110	14 x 100
Tubes supplied	Yes	Yes	Yes	Yes	Yes

Swing Out Rotor R-235 M (8 x 96 WELL)		
C-23 Plus	CPR-24i Plus	C-30 Plus CPR-30 Plus
 Max. Speed RPM 3200 Max RCF (g) 1370		
Image of Tube		
Tube Capacity (ml)	8 x 96 wells	
Type of Tube	Microtiter	
No of tubes/rotor	8	
W x D x L (mm)	127 x 84 x 18	
Tubes supplied	No	

Swing Out Rotor R-230 BM + R-233 BM (108 x 5ml (RIA))	
C-23 Plus	C-30 Plus CPR-30 Plus
Max. Speed RPM 5100 Max RCF (g) 3980	Max. Speed RPM 4500 Max RCF (g) 3500
	
Image of Tube	
Tube Capacity (ml)	108 X 5 (RIA)
Type of Tube	RIA Plastic Tubes
Insert or Cushion No.	No
No of tubes/Carrier	27
No of tubes/rotor	108
Ø x L (mm)	12.4 x 76
Tubes supplied	Yes

Swing Out Rotor R-240 M (16 x 50ml)	
C-30 Plus CPR-30 Plus	
Max. Speed RPM 3200 Max RCF (g) 2110	
	
Image of Tube	
Tube Capacity (ml)	50
Type of Tube	Conical
Insert or Cushion No.	C 50F Included
No of tubes/Insert	4
No of tubes/rotor	16
Ø x L (mm)	29 x 116
Tubes supplied	No

Swing Out Rotor R-240 AM (28 x 15ml)	
C-30 Plus CPR-30 Plus	
Max. Speed RPM 3200 Max RCF (g) 2110	
	
Image of Tube	
Tube Capacity (ml)	15
Type of Tube	Conical
Insert or Cushion No.	C 15F Included
No of tubes/Insert	7
No of tubes/rotor	28
Ø x L (mm)	16.5 x 120
Tubes supplied	No

Swing Out Rotor R-246 AM (4 x 25ml)			Image of Tube	
C-24 Plus CPR-24 Plus	CPR-24i Plus	C-30 Plus CPR-30 Plus	Tube Capacity (ml)	25
 <p>Max. Speed RPM 11000 Max RCF (g) 12180</p>			Type of Tube	Glass
			Insert or Cushion No.	C 25R Included
			No of tubes/Insert	1
			No of tubes/rotor	4
			Ø x L (mm)	28 x 62.5
			Tubes supplied	Yes