



Homogenizers

RQ - 127A / D suitable for homogenizing of tissues at high speed available. Chuck to hold shaft upto 10 mm dia. Fitted with PTFE grinding pestle. Glass homogenising cup with rod of 50ml supplied as standard, optional available 5 ml, 15 ml, 25 ml & 100 ml. RQT-127A/D is with digital RPM indicator

RQ-127 A / D

Technical Data								
Model	Stirring Shaft		Impeller	Motor	Speed (rpm)		Stirring Capacity	Digital Speed Indicator
	Length (mm)	Ø (mm)	Type		Min	Max	Max Volume	
RQ-127A/D	275	8	Teflon Pestle	AC/DC 1/8 HP	300	8000	upto 100 ml.	✕
RQT-127A/D	275	8	Teflon Pestle	AC/DC 1/8 HP	300	8000	upto 100 ml.	LED

Supply: 220-240 Volts 50 Hz Single Phase.