

Salient Features:

- Designed for mixing liquids in Schools, Laboratories & Factories
- Touch/Continuous Operation mode Selection through bi-directional Switch
- Speed Regulation through knob
- Interchangeable mixing heads for use with variety of tubes
- Supplied with all interchangeable mixing heads



CM-101 Plus



CM-101

Technical Data		
Model	CM-101 Plus	CM-101
Shaking Movement	Orbital	
Orbital Diameter (mm)	4	
Motor Type	Shaded-Pole Motor	
Permissible ON time	100% Power 30mins	
Speed range (rpm)	Upto 2900	
Run type	Continuous / touch operation	
Dimensions (WxDxH) mm	185 x 150 x 165	178 x 141 x 148
Weight (Kg)	2.8	3.6
Supply: 220-240 Volts 50 Hz Single Phase.		

Interchangeable Mixing Heads	
Model	Description
CM-01	Universal top 100 mm- plate attachment
CM-01/01	Tube adapter for 18 holes test tubes Ø 10mm [use with CM-01]
CM-01/02	Platform pad for < Ø 99mm tubes and small vessels [use with CM-01]
CM-01/03	Tube adapter for 8 holes test tube Ø 20mm [use with CM-01]



CM - 01



CM - 01/01



CM - 01/03