




Salient Features for CM-101 Ultra

- Compact size with high capacity
- Stable operation with heavy base
- Brushless motor with strong power
- Multipurpose adaptor / large handling capacity
- Feather touch keys
- Easy to switch, touch/run

CM-101 Ultra

 Hybrid Adaptor Max speed RPM 3000				
Image of Tube				
Tube Capacity (ml)	2 / 1.5	0.5	2 Strips of 8 x 0.2	1 x 96 wells
Type of Tube	Micro	Micro	PCR	MTP Plate
No of tubes	8	12	2	1
Ø x L (mm)	11 x 41	7.5 x 31	6 x 78	W127 x D84 x L18
Tubes supplied	No	No	No	No

Technical Data	
Model	CM-101 Ultra
Speed Range	100-3000 rpm
Control Mode	Continuous / Touch
Motor	Brushless motor
Amplitude	4.5mm (Orbital)
Display of Speed	LCD
Timer range	1 min - 8 hours
Shell Protection Level	IP-21
Dimension (WxDxH)	150x150x77mm

Supply: 220-240 Volts, 50 Hz Single Phase.