





RM-03 Pro

Salient Features for RM-03 Pro:

- 9 Memory programs
- Safety lid lock to prevent running unless lid is closed
- LCD display of parameters
- BLDC Motor
- Supplied with 2 Rotors
- Safety over temperature cut-off for sample protection

Rotors for RM-03 Pro				
 Micro Tube Rotor Max speed RPM 14500 Max RCF (g) 14104		 4 Strips PCR Rotor Max. Speed RPM 14500 Max RCF (g) 12576		
Image of Tube				Image of Tube
Tube Capacity (ml)	2	0.5	0.2	Tube Capacity (ml)
Type of Tube	Micro	Micro	Micro	Type of Tube
Insert Supplied	No	Yes	Yes	Type of Tube
No of tubes/Insert	-	1	1	No of tubes/rotor
No of tubes/rotor	12	12	12	Ø x L (mm)
Ø x L (mm)	11 x 41	7.5 x 31	6 x 22	Tubes supplied
Tubes supplied	No	No	No	No

(Rotor supplied with all suitable inserts)

Technical Data	
Model	RM-03 Pro
Max Speed	14500rpm
Max. RCF	14100g
Rotor Supplied	Micro tube rotor, PCR strip rotor
Timer range	30secs-59min
No. of Programs	9
Dimension (WxDxH) (mm)	180x270x105

Supply: 220-240 Volts, 50 Hz Single Phase.



RM-03 Plus

Salient Features:

- Microprocessor controlled with last run memory feature
- Large back - lit LCD display shows multiple

parameters at a time

- Brushless DC motor - maintenance free
- Imbalance detection with auto cut off
- Lid lock safety
- Intuitive simple interface for quick and convenient setting
- 1 to 999 Mins countdown timer with Min sec display
- One touch Short- Spin operation
- Unique airflow design - ensures lesser heat generation and lesser noise
- Rotor, inner centrifuge lid and lid latch are made of metal for maximum operational safety
- Soft - touch lid closure for ergonomic lid locking
- Settable RPM/RCF

Technical Data	
Specifications	RM-03 Plus
Rotor Capacity	12 place closed rotor with metal safety lid
Maximum Speed RPM / RCF (g)	15000 rpm / 15596 g
Imbalance detection	Imbalance Detection with auto cut - off
Speed / Time setting	500 to 15000 RPM / Min-Sec timer
Dimensions (WxDxH) in mm	230 x 262 x 131 mm
Rotors and Adaptors (included in standard pack)	12 place x 1.5/2.0 ml microtubes
	Reduction Adaptors for 0.4/0.5 ml microtubes
	PCR Strip Rotor (2 x 8 x 0.2 ml with metal safety lid)
Supply: 220-240 Volts, 50 Hz Single Phase.	