



- **Salient Features for RM-02 Pro**
- Higher Speed upto 7000 RPM
- Multi-rotor, more working capacity
- Low Centre of gravity for stable operation
- High frequency and wide voltage
- Low noise
- Supplied with Hybrid Rotor-accomodates microtubes & PCR Strips



RM-02 Pro

Rotor for RM-02 Pro				
		Micro Hybrid Rotor Max speed RPM 7000 Max RCF (g) 3280		
Image of Tube				
Tube Capacity (ml)	2 / 1.5	0.5	0.2	4 Strips of 8 x 0.2
Type of Tube	Micro	Micro	Micro	PCR
Insert Supplied	No	Yes	Yes	No
No. of tubes/Insert	-	1	1	-
No. of tubes/rotor	12	12	12	4
Ø x L (mm)	11 x 41	7.5 x 31	6 x 22	6 x 78
Tube supplied	No	No	No	No

(Rotor supplied with all suitable inserts)

Technical Data	
Model	RM-02 Pro
Control Mode	Digital Display & Button
Max Speed	7000rpm
Max RCF	3280xg
Timer Range	10sec-99min
Dimension (W x D x H) (mm)	180 x 185 x 135

Supply: 220-240 Volts, 50 Hz Single Phase



RM-02 Plus

Powerful and easy-to-use **RM 02 Plus mini centrifuges** are small enough so that each workstation can be equipped with a personal centrifuge. Extremely compact footprint saves valuable bench space

Salient Features :

- Brushless DC (BLDC) Motor
- Settable speed (upto 6000 RPM) and time (upto 99 minutes). digital display can process upto 8 microtubes of 1.5/2 ml at a time
- Replaceable strip rotor & adaptors for 0.2 ml & 0.4/0.5 ml microtubes

Technical Data	
Specifications	RM-02 Plus
Rotor Capacity	8 place closed rotor for 1.5/2.0ml microtubes
Speed RPM / RCF (g)	6000/2000 g
Digital Display	Yes
Dimensions (WxDxH) in mm	162 x 157 x 115 mm
Rubber suction pads	Yes
Rotors and Adaptors (included in standard pack)	8 Slot x 1.5/2.0 ml Closed Rotor PCR Strip Rotor (2 x 8 x 0.2 ml) Reduction Adaptors for 0.2 / 0.4 ml microtubes
Supply: 220-240 Volts, 50 Hz Single Phase.	