



Mini Magnetic Stirrers

Salient Features:

- LED Display of Speed
- Compact unit, small foot print
- Last memory recall
- Small wall mounted power adaptor
- IP-42 Protection level
- Supplied with PTFE Stirring Paddle



MS-1 Plus

Technical Data	
Model	MS-1 Plus
Max Stirring Volume (Water)	1 Litre
Motor Type	Brushless
Speed range (RPM)	200-2000
Working Plate dia (mm)	135
Stirring paddle (PTFE coated) (mm)	11 x 30
Shell Protection Level	IP-42
Dimension (WxDxH) (mm)	168x160x47

Supply: 220-240 Volts, 50 Hz Single Phase.

Specifications subject to change without notice

