



RS-10 Pro

Salient Feature:

- Timing Function
- Adjustable speed
- Speed increases in step of 1 RPM
- Choice of timed & Continuous mode
- Suitable for a variety of containers

Technical Data	
Model	RS - 10 Pro
Motor	Geared
Speed	10 - 80rpm
Timer Range	1min - 24hrs
Max Load	2kg
Display	LCD
Tilt Angle	7°
Tray size in mm	Silicone tray 250 x 250
Dimensions (WxDxH)	355 x 260 x 135 mm

Supply: 220-240 Volts, 50 Hz Single Phase