

Rotor Listing for R-8C

Models	Capacity	Max. Speed RPM	Max. RCF 'g'	Max. Speed RPM	Max. RCF 'g'	Models	Capacity	Max. Speed RPM	Max. RCF 'g'	Max. Speed RPM	Max. RCF 'g'
		R-8C BL		R-8C				R-8C BL		R-8C	
Swing out Rotors (with graduated glass tubes)}						Angle Rotors (with polypropylene tubes)}					
R-81	16x15 ml	4500	3485	4000	2750	R-83	16x15 ml	6000	5000	5250	3800
R-81 A	12x15 ml	4500	3485	4000	2750	R-83 A	12x15 ml	6000	5000	5250	3800
R-81 B	8x15 ml	4500	3440	4200	3000	R-83 B	8x15 ml	6000	4700	5250	3600
R-82	6x50 ml	4500	3530	4000	2800	R-84	6x50 ml	6000	4870	5250	3730
R-82 A	4x50 ml	4500	3585	4000	2750	R-84 A	4x50 ml	6000	4730	5250	3620
R-89 (With Biosafety Cap)	4x100 ml	4000	2500	-		R-89	24x15 ml	6000	5070	-	-
						R-88	4x100 ml	6000	4870	Discontinued	

Reduction Adaptors for R-8C

Models	Capacity	R-81	R-81A	R-81B	R-82	R-82A	R-83	R-83A	R-83B	R-84	R-84A	R-88	R-89	R-90
RA-151	1 x wintrob	Discontinued												
RE 15-3R (OLD - RA-152)	1 x 3ml	✓	✓	✓	✗	✗	✓	✓	✓	✗	✗	✗	✓	✗
RE 15-5R (OLD - RA-153)	1 x 5ml	✓	✓	✓	✗	✗	✓	✓	✓	✗	✗	✗	✓	✗
RE 15-8R (OLD - RA-154)	1 x 8ml	✓	✓	✓	✗	✗	✓	✓	✓	✗	✗	✗	✓	✗
RA-155	1 x 4ml	Discontinued												
RE 15-6R (OLD - RA-156)	1 x 6ml	✓	✓	✓	✗	✗	✓	✓	✓	✗	✗	✗	✓	✗
RE 100-50 (OLD - RA-157)	1 x 50ml	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓
RA-158	1 x 30ml	Discontinued												
RA-160	1 x 15ml	Discontinued												
RE 100-50RR (OLD - RA-161)	1 x 50ml	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✗	✗
RE 100-15R(3) (OLD - RA-166)	3 x 15ml	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✗	✓
RE 50-15RR (OLD - RA-501)	1 x 15ml	✗	✗	✗	✓	✓	✗	✗	✗	✓	✓	✗	✗	✗
RA-502	1 x 30ml	Discontinued												
Suitable for Falcon Tubes (Without Tubes)														
RE 50-15FR (OLD - RA-501F)	1 x 15ml	✗	✗	✗	✓	✓	✗	✗	✗	✓	✓	✗	✗	✗
RE 100-15FR (OLD - RA-160F)	1 x 15ml	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✗	✓
RE 100-50FR (OLD - RA-161F)	1 x 50ml	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✗	✓
Suitable for Vacutainer Tubes (Without Tubes)														
RE 15-5/7SR (OLD - RA-162)	1 X upto 5ml (75 mm long)	✓	✓	✓	✗	✗	✗	✓	✗	✗	✗	✗	✗	✗
RE 15-5/7LR (OLD - RA-163)	1 X upto 6ml (100 mm long)	✓	✓	✓	✗	✗	✗	✓	✗	✗	✗	✗	✗	✗
RE 50-10R (OLD - RA-164)	1 X upto 10ml (100 mm long)	✗	✗	✗	✓	✓	✗	✗	✗	✓	✓	✗	✗	✗
RE 100-5/7L(4) (OLD - RA-167)	4 X upto 6ml (100 mm long)	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✗	✓
RE 100-10(3) (OLD - RA-168)	3 X upto 10ml (100 mm long)	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✓	✗	✓

Rotors Listing for RM-12C

Models	Capacity	Max. Speed RPM	Max. RCF 'g'
RM-1210	Micro Hematocrit Rotor for 24 Capillaries (75 X 1 mm)	12000	15300
RM-1211	Micro Hematocrit Rotor for 12 capillaries (QBC -II Plus)	12000	15300
RM-1212	Heparinized capillary tubes (75mm x 1 mm bore)		
RM-1213	Reading Device for RM-1210		
Angle Rotors			
RM-1214	10 X 2ml with Glass Tubes & reduction adaptors of 1ml & 0.5ml	16000	16600
RM-1215	24 X 1.5ml with Tapered bottom Polypropylene tubes & reduction adaptors of 0.5ml & 0.2ml	11850	12715
RM-1216	8 X 5ml with round bottom Polypropylene tubes & reduction adaptors of 4ml & 2ml	11850	13500
RM-1217	4 strips of 8 X 0.2ml (PCR Strip Tubes)	12000	11910

Specifications subject to change without notice