



Smart Vortex Mixers

Salient Features for CM-101 Pro:

- Automatic start / stop on sensing operator hand
- Soft silicon cap
- Visual indication of mixing
- Large mixing area
- Feather touch keys
- Adjustable Speed
- Suitable for a variety of containers



CM-101 Pro

Technical Data	
Model	CM-101 Pro
Speed Range	100-3000 rpm
RPM Increment	100 rpm
Control Mode	Continuous / Auto sensing
Motor	Brushless motor
Amplitude	4mm (Orbital)
Display of Speed	LCD
Dimension (WxDxH)	150x150x45mm

Supply: 220-240 Volts, 50 Hz Single Phase.

Specifications subject to change without notice

