

**RSB - 12****Water Bath Shaker RSB - 12**

These Water bath Shakers have wide applications in Biotechnology, Micro - Biology, Chemistry, Quality Control and other life sciences laboratories.

The bath has size of 355 mm (W) x 300 mm (D) x 165 mm (H). A special pyramid type cover provided avoids evaporation of water and makes condensed water to fall back in the bath. A cock is provided to facilitate draining of Bath contents

The removable stainless steel tray houses lotus clamps for variety of flask capacities. **Standard unit supplied with tray accommodating 12 flasks x 100 ml.**

Following capacity are available optionally.

50 ml x 12 flasks	150 ml x 9 flasks
250 ml x 6 flasks	500 ml x 4 flasks
1000 ml x 2 flasks	2000 ml x 1 flask
Test Tube racks	

Salient Features:

- Inner chamber and outer body of SS 304
- Heating up to 95°C
- Digital Temperature Indicator cum Controller with precise accuracy of $\pm 1^{\circ}\text{C}$
- Shaking speed range from 40 to 180 strokes per min
- Stroke length 25 mm

Supply: 220-240 Volts 50 Hz Single Phase