



Water Baths



RWB-6



RWB-12

REMI Water Bath is multipurpose water bath with six holes/ twelve holes to hold different size samples for easy & hygienic heating. Water Bath provides precise and constant temperature helps achieving uniform heating for the samples.

Salient Features:

- Microprocessor control with timing function
- Digital display
- Audible and visual alarm for over temperature
- **GMP versions with Inner & Outer SS 304**

Technical Data		
Parameter	RWB - 6	RWB - 12
Temperature Range	RT $\pm$ 5°C ~99°C	RT $\pm$ 5°C ~99°C
Temperature Stability	$\pm$ 0.5°C	$\pm$ 0.5°C
Deviation Ability	$\pm$ 2°C	$\pm$ 2°C
Chamber Volume	15 Ltrs	28 Ltrs
Interior Dimension (W x D x H) (mm)	450 x 300 x 110	500 x 390 x 125
Timing Range	1~999 min.	1~999 min.
No. of hole	Six holes	Twelve holes
Heater Capacity (Watts)	1500	3000
Supply: 220-240 Volts, 50 Hz Single Phase.		

Specifications subject to change without notice

