

**ZERO
COMPROMISE MODELS**

- BRUSHLESS
- ECONOMICAL



RQ-5 Plus



RQ-100 Plus

Direct Drive Stirrers

REMI next generation of laboratory stirrers for stirring / mixing etc of chemical, pharmaceutical, food and cosmetic products.

Salient Features :

- High torque even at low speeds
- Monitoring of set & actual speed
- Wide range of impellers
- Advanced Stirrers with microprocessor technology
- Advanced Stirring speed control
- Constant speed even with change of load / Voltage
- Easy -to- read display of set & actual parameters
- Last parameter recall. Ideal for repetitive processes
- Motor protection in case of continuous overload
- Special base design for topple free operation
- Stirring Shaft SS316, Retort Rods SS304

Technical Data

Features	RQ-5 Plus	RQ-100 Plus
Max. stirring capacity (Water) (Liters)	5	100
Motor Type	BLDC	AC Induction, 3 phase
Motor Rating	10 watt	100 watt
Motor Range (rpm)	100 -1500	50-4000
Speed Display	LED	LED
Timer	0 to 99 Minutes	0 to 99 Minutes
Motor Overload Protection	LED light flash, auto stop	LED light flash, auto stop
Chuck Ø (mm)	1.5 to 10	1.5 to 13
Stirring Shaft Ø x length (mm)	6 x 380 mm	12 x 500
Impeller Type / Ø (mm)	Axial Turbine / 36	Pitched x 2 / 65
Dimensions [W x D x H] (mm)	80 x 180 x 370	320 x 285 x 620
Stand length (mm)	500	1200
Weight (Kg)	6.5	19
Power (W)	10 watt	175 watt
Supply: 220-240 Volts, 50 Hz Single Phase.		

Specifications subject to change without notice

