



RCB-13

Provided for precise and constant temperature and auxiliary heating in colleges, industrial and mining enterprises and scientific research departments.

- **Salient Features:**
- LCD Microprocessor temperature controller
- Timing function Range 99.59 hrs
- Rounded corner SS 304 chamber
- Audible and visible alarm for temperature and water level
- With interface to connect externally
- Drain valve for faster discharge & easy cleaning

Technical Data	
Model	RCB - 13
Temperature Range	-20~100°C
Precision	± 0.2°C
Chamber Volume	13L
Pump (flow rate)	8 Litre/min
Chamber Dimension (mm)	240 x 170 x 200
External Dimensions (mm)	380 x 600 x 875

Supply: 220-240 Volts, 50Hz Single Phase.

Specifications subject to change without notice