

**REMI**

# Dry Bath Shakers

**NEW**  
INTRODUCTION



**DBS-10**



**Heat Lid**

**REMI DRY BATH SHAKING INCUBATOR** adopts advanced design concept and manufacturing technology. It has the advantages of compact and firm structure, high versatility and simple operation. **DRY BATH SHAKING INCUBATOR** has multi-purpose functions such as cooling, heating and shaking. It is ideal for sample incubation, catalysis, mixing and storage. A wide variety of **replaceable aluminium blocks** (to be ordered separately) for different types of tubes make the unit truly versatile. **Cooling is achieved by Peltier technology.** Audio indication on completion of cycle. **Safety protection for over speed running by Auto detection of block height.**

## Technical Data

| Model                                                  | DBS-10                                                         |
|--------------------------------------------------------|----------------------------------------------------------------|
| Temperature Range                                      | Room temp. minus 15°C to 100°C                                 |
| Temperature Stability                                  | ±0.5°C                                                         |
| Max heating rate                                       | 5°C / min                                                      |
| No. of Programs                                        | 9                                                              |
| Sequential Program Run                                 | Yes                                                            |
| Speed                                                  | 100-1500 rpm(Low Block) / 100-800 rpm(High Block)              |
| Timer Range                                            | 5secs - 24hours                                                |
| Heat Lid for improved Temperature stability (Optional) | Suitable for MS-96x0.2, MS-24x0.5, MS-24x1.5, MS-24x2, MS-Flat |
| Dimension (WxDxH)                                      | 190x265x150mm                                                  |

Supply: 220-240 Volts, 50Hz Single Phase.

## Essential Blocks for DBS-10 (To be ordered separately)

| Block Model No. | Block Height | Capacity                                 |
|-----------------|--------------|------------------------------------------|
| MS-96x0.2       | Low          | 96 x 0.2ml tubes/96-well plate/PCR Strip |
| MS-24x0.5       | Low          | 24 x 0.5ml tubes                         |
| MS-24x1.5       | Low          | 24 x 1.5ml conical bottom tubes          |
| MS-24x2         | Low          | 24 x 2ml round bottom tubes              |
| MS-8x15         | High         | 8 x 15ml conical tubes                   |
| MS-4x50         | High         | 4 x 50ml conical bottom tubes            |
| MS-Flat         | Low          | Suitable for 96 well Elisa/Culture plate |



MS-96x0.2



MS-24x0.5



MS-24x1.5




MS-24x2



MS-8x15



MS-4x50



MS-Flat

Specifications subject to change without notice

