

**WORLD CLASS
PREMIUM RANGE**



RS - 24 Plus

Ideal for mixing and development of cultures, chemicals, solvents, and assays etc. in Microbiological, Cell Culture & Life Science laboratories.

Brushless Induction drive motor with frequency drive makes the unit suitable for continuous non-stop operation. Step less electronic frequency control ensures gentle start and maintains preset speed. Compact counter balanced drive mechanism ensures high stability and reliability even in continuous operation & uneven load distribution.

Salient Features :

- 4" Large LCD display
- Universal Platform with rounded corners to accommodate interchangeable clamps of assorted sizes
- Brushless induction motor with VFD speed control
- Elegant Plastic front
- Variable speed 20 -300 RPM
- Shaking amplitude 40 mm
- Auto start after power resumption

Technical Data

Models	Platform Size	Max. Volume	Max. No. of clamps (To be ordered separately)					
			2000 ml.	1000 ml.	500 ml.	250 ml.	150 ml.	100ml/ 50ml
RS-24 Plus	24" x 24"	26 Litres	13	16	25	36	64	64

Supply: 220-240 Volts, 50 Hz Single Phase.

Specifications subject to change without notice