



The REMI Phlebotomy Chairs



Flabo X



Flabo XR

REMI **Flabo X** and **Flabo XR** are advanced phlebotomy chairs designed for efficient sample collection and patient comfort.

While **Flabo X** is a fixed-position chair offering sturdy support and ergonomic armrest design, **Flabo XR** is a high-end recliner chair with independently adjustable leg and back sections through a dual hydraulic system.



Additional Features of Flabo XR

- Equipped with a dual hydraulic system for independent leg and back section movement using separate levers.
- Reclining function allows transition from sitting to flat bed position, supporting patient assessment and comfort.
- Mounted on sturdy dual-lock PU castor wheels for emergency mobility.

Salient Features

Comfort

- Ergonomic design to ensure patient comfort during phlebotomy.
- Broad sitting space with non-deformable PU foam accommodates even obese patients (up to 150 kg).

Support

- Broad, firm, and foamed handrest for enhanced arm support during blood collection.
- Handrest are adjustable vertically, horizontally, forward and backward to help achieve ideal phlebotomy positioning.

- Soft cushioned headrest adds extra support for reclining patients.
- Robust steel frame construction ensures long-term durability and stability in both models.

Hygiene

- Stain-resistant and washable PU upholstery.
- Smooth and rounded edges make cleaning simple and efficient.

Technical Specifications

Parameters	Flabo XR	Flabo X
Overall Dimensions (L X W X H)	1791 x 533 x 516 mm	760 x 770 x 1230 mm
Weight	35 kg	25 kg
Possible postures	Flat bed & Sitting posture	Sitting posture
Hydraulic system	2 (for Leg & Back section)	-
Load caring capacity	150 kg (max)	150 kg (max)