



NEYA 2

Simplifies Day to Day Work





Image of Swing Out Rotors					
Model	SM 4-50	SC 8-15	SM 6-15	S 16-5	
Max speed RPM	4500	4500	4500	4500	
Max RCF (g)	2810	2880	3000	2950	
Image of Tube					
Tube Capacity (ml)	50	15	15	5/7	5/7
Type of Tube	Glass	Glass	Conical	Blood Collection Tube (Long)	Blood Collection Tube (Short)
Insert or Cushion No.	C 50R Included	C 15R Included	C 15F	Included	Included
No of tubes/Insert	1	1	1	1	1
No of tubes/rotor	4	8	6	16	16
Ø x L (mm)	30 x 110	17 x 110	16.5 x 120	12.5 x 106	12.5 x 81
Tubes supplied	Yes	Yes	No	No	No

Image of Angle Rotors					
Model	AM 6-50			AM 12-15	AM 8-15
Max speed RPM	4500			4500	4500
Max RCF (g)	2740			2760	2760
Image of Tube					
Tube Capacity (ml)	50	50	15	15	15
Type of Tube	Plastic	Conical	Conical	Glass	Conical
Insert or Cushion No.	C 50R Included	C 50F	RE 50-15FR (OLD - RA-501F)	C 15R Included	C 15F
No of tubes/Insert	1	1	1	1	1
No of tubes/rotor	6	6	6	12	8
Ø x L (mm)	30 x 100	29 x 116	16.5 x 120	17 x 110	16.5 x 120
Tubes supplied	Yes	No	No	Yes	No