



NEYA 2

Simplifies Day to Day Work



**Specialized Rotor
for Diagnostic labs**



S 16-5

Specially Designed Hybrid Rotor
(Long & Short Blood Tubes -16x 5/7ml)

REMI Neya 2 with European Design is a perfect combination of performance & Characteristics. Neya 2 is available with brush less motor, door safety interlock, digital speed indicator and timer. Ideal for routine centrifuging requirements of medical practitioners, diagnostic laboratories and hospitals. Precise control on speed is through microprocessor based VFD controller, helps in efficient results for sensitive diagnostic tests.

REMI Neya 2 is designed with compact footprint considering the day to day requirements of small to medium medical diagnostic laboratories.

Specialized hybrid rotor to fit both popular sample tubes simultaneously.

1. Short Blood collection Tubes (Size - 12.5mm X 81 mm)
2. Long Blood collection Tubes (Size - 12.5mm X 106 mm)

Technical Specifications		
Particulars		Neya 2
Capacity	ml	300 ml
Max. Speed	RPM	4500
Max. RCF	'g'	3000
W X D X H	mm	320 x 420 x 280

Supply: 220 V, 50 Hz Single phase





NEYA 2

Simplifies Day to Day Work



Salient Features

- Brushless induction motor
- Speed control through VFD
- Microprocessor controller
- Feather touch keys for Smooth and soft start
- Green LED display with last parameter recall
- Maximum RPM 4500 and maximum RCF 3000 g
- Lid interlock, Dynamic brake and Imbalance detection
- Emergency lid opening
- Accurate speed ± 10 RPM
- Timer up to 99 minutes
- Calibration window
- Choice of swing-out & angle rotors with special blood tube rotor






Image of Swing Out Rotors				
Model	SM 4-50	SC 8-15	SM 6-15	S 16-5
Max speed RPM	4500	4500	4500	4500
Max RCF (g)	2810	2880	3000	2950
Image of Tube				
Tube Capacity (ml)	50	15	15	5/7
Type of Tube	Glass	Glass	Conical	Blood Collection Tube (Long)
Insert or Cushion No.	C 50R Included	C 15R Included	C 15F	Included
No of tubes/Insert	1	1	1	1
No of tubes/rotor	4	8	6	16
Ø x L (mm)	30 x 110	17 x 110	16.5 x 120	12.5 x 106
Tubes supplied	Yes	Yes	No	No

Image of Angle Rotors			
Model	AM 6-50	AM 12-15	AM 8-15
Max speed RPM	4500	4500	4500
Max RCF (g)	2740	2760	2760
Image of Tube			
Tube Capacity (ml)	50	15	15
Type of Tube	Plastic	Conical	Glass
Insert or Cushion No.	C 50R Included	C 50F	C 15R Included
No of tubes/Insert	1	1	1
No of tubes/rotor	6	6	12
Ø x L (mm)	30 x 100	29 x 116	17 x 110
Tubes supplied	Yes	No	Yes