


Blood Donor Chair



- Contoured for better spine support, increases donor comfort & safety
- Natural feet high position increases blood flow, supports donation process
- High grade modulated FRP body (Fibre Reinforced Polymer) on steel frame
- Motorized remote slide control for head low - feet high (vice versa) position
- Non deformable PU foam with lumbar support grooves
- Long lasting eco-leather upholstery for easy maintenance
- Broad & adjustable hand rests supports phlebotomy
- Lockable castors to ensure firm ground holding
- Both side adjustable IV stand with 3 trays for accessories
- Interface for blood collection monitor

Technical Specifications

Parameters	Dzire 1X
Seat width	652 mm
Length	1500 mm
Height max head side	1210 mm
Height max leg side	1050 mm
Seat height from floor	635 mm
Weight	75 kg
Weight capacity	200 kg
Hand rest dimension (W x L mm)	560 x 235
Compliance	IEC 60601-1, IEC 60601-1-2

Supply : 220-240 volts, 50 Hz, SINGLE PHASE

