



Mobile Transport Box MBT Series



MBT 100

MBT 100 is a refrigerated mobile transport box for easy & safe transfer of temperature sensitive biological products.

It provides active temperature control with a compact powerful compressor. The flexibility of operations with options of AC / DC power makes it an ideal choice in maintaining cold chain on the move.

The wide temperature setting range supports requirement of different products. The display & alarm system control helps in continuous monitoring of temperature.

Salient Feature

- Digital display of temperature
- High / low temperature alarm
- Wide temperature range (-20°C to +24°C)
- CFC Free Refrigerant R134a
- Lockable castor wheels with recessed handle
- SS Hinges & pad lockable door latch for easy closure of lid

Technical Specifications

Parameters	Unit	MBT 100
Dimension (L x W x H)	mm	500 x 825 x 635
Basket Shelves	3 Nos	One small & two regular baskets
Capacity	Ltr	100
Weight	Kg	30
Temperature range	°C	-20°C to + 24°C
Voltage	V	DC : 12 / 24 V & AC : 220-240V, 50Hz
Body Construct	Powder coated CRCA Steel exterior & AL inner lining	