



## Portable Donor Chair



**PDC 1X**



**Folded**

REMI Portable Donor Chair model PDC 1X, specially designed for blood donation camps to provide comfort and safety to the blood donors.

- Ergonomic robust design
- Stable & aesthetically beautiful
- Complete construct with stainless steel
- High quality anti bacterial fabric
- Foldable for easy transport
- Adjustable safety belt for patient safety
- Towel bar at top
- 30 kg light weight for easy transportation
- IV Stand optional (standalone)
- Comfortable seating
- Trendelenburg (Head low) & Reverse-Trendelenburg (Sitting Posture) position
- Elegant ergonomic knobs for the hand rest adjustment
- Adjustable head rest positions for easy phlebotomy
- Elegant ergonomic knob for hand rest adjustment
- Adjustable hand rest positions (Up / down and horizontal / circular)
- Comfortable molded foamed hand rests
- Fixed Castor Wheels attached for effortless handling and easy mobility

### Technical Specifications

| Parameters                     | PDC 1X    |
|--------------------------------|-----------|
| Seat width                     | 500 mm    |
| Length                         | 1740 mm   |
| Seat Height from floor         | 1270 mm   |
| Weight Capacity                | 150 kg    |
| Hand rest dimension (W x L mm) | 150 x 450 |
| Net Weight                     | 30 kg     |