



Automated Plasma Expresser



XPRESS - A

XPRESS - A is a semi-automated plasma expressor for quick & easy component separation. The optical sensing of RBC's & automatic clamping reduces manual error & prevents RBC contamination.

Salient Features

- Automated fast acting clamp mechanism, no RBC contamination
- Spring loaded transparent compression plate for quick plasma expression
- Adjustable spring tension
- Ergonomic handle for easy operation
- LED power indication, buzzer alarm on clamping
- Suction holders for firm base grip
- Flexibility of manual operation



Control Panel

Technical Specifications

Parameters	Unit	Xpress A
Dimensions (W x D x H)	mm	198 x 245 x 309
Max Tube Diameter	mm	6.5
Weight	kg	4
Operating Temperature	°C	0 to 40
Input Supply		230 V, 50 Hz