



Manual Plasma Expressor



XPRESS - M

XPRESS - M is specially designed plasma expressor for separating out the blood components. The front panel is loaded with high tension spring for quick transfer of plasma into other bag.

XPRESS - M is ergonomically designed for long lasting trouble free operations

Salient Features

- User friendly with see through front panel to allow visibility of layer movement
- Spring loaded transparent compression plate for quick plasma expression
- Adjustable spring tension
- Ergonomic handle for easy operation
- Suction holders for firm base grip
- Metallic base covered with rounded ABS body
- Rust free, easily cleanable surface

Technical Specifications

Parameters	Unit	Xpress M
Dimensions (W x D x H)	mm	175 x 230 x 275
Weight	kg	3.5
Body Construct		Powder Coated CRCA Steel
Front Panel		Acrylic