




Angle Rotors for R-12C Plus

Angle Rotor AR 12-2 Max speed RPM 16000 Max RCF (g) 17200 	Image of Tube				
	Tube Capacity (ml)	2	1	0.5	0.2
	Type of Tube	Micro	Micro	Micro	Micro
	Insert or Cushion No.	No	RE 2-1	RE 2-05	RE 2-02
	No of tubes/Insert	-	1	1	1
	No of tubes/rotor	12	12	12	12
	Ø x L (mm)	11 x 41	7.5 x 49	7.5 x 31	6 x 22
	Tubes supplied	No	No	No	No


(Rotor supplied with all suitable inserts)

Angle Rotor AR 24-2 Max speed RPM 15000 Max RCF (g) 21000 	Image of Tube			
	Tube Capacity (ml)	2	0.5	0.2
	Type of Tube	Micro	Micro	Micro
	Insert or Cushion No.	No	RE 2-05	RE 2-02
	No of tubes/Insert	-	1	1
	No of tubes/rotor	24	24	24
	Ø x L (mm)	11 x 41	7.5 x 31	6 x 22
	Tubes supplied	No	No	No




(Rotor supplied with all suitable inserts)


Angle Rotor AR 28-1.5 Max speed RPM 16000 Max RCF (g) 25900 	Image of Tube			
	Tube Capacity (ml)	1.5	0.5	0.2
	Type of Tube	Micro	Micro	Micro
	Insert or Cushion No.	No	RE 2-05	RE 2-02
	No of tubes/Insert	-	1	1
	No of tubes/rotor	28	28	28
	Ø x L (mm)	11 x 41	7.5 x 31	6 x 22
	Tubes supplied	No	No	No

(Rotor supplied with all suitable inserts)


Angle Rotor AR 8-5 Max speed RPM 11850 Max RCF (g) 13500 	Image of Tube			
	Tube Capacity (ml)	5	4	2
	Type of Tube	Plastic	Plastic	Micro
	Insert or Cushion No.	No	RE 5-4	RE 5-2
	No of tubes/Insert	-	1	1
	No of tubes/rotor	8	8	8
	Ø x L (mm)	12 x 76	10 x 75	11 x 41
	Tubes supplied	No	No	No

(Rotor supplied with all suitable inserts)

Image of Angle Rotors			
Model	AR 24-75 HE	AR 12-75 HE	P 4-8
Max speed RPM	12000	12000	15000
Max RCF (g)	15300	15300	21000
Image of Tube			
Tube Capacity (ml)	-	-	4 strips of 8 PCR
Type of Tube	Capillary	Capillary QBC-II Plus	PCR
No of tubes/rotor	24	12	4
Ø x L (mm)	1 x 75	-	6 x 78
Tubes supplied	Yes	No	No

Angle Rotor AR 12-5 Max speed RPM 13000 Max RCF (g) 16000 	Image of Tube	
	Tube Capacity (ml)	5
	Type of Tube	Micro
	Insert or Cushion No.	No
	No of tubes/Insert	-
	No of tubes/rotor	12
	Ø x L (mm)	16 x 59
	Tubes supplied	No

Reading Device RD 24-75



W x D x H

127 x 84 x 18

Specifications subject to change without notice