



Medico Doctor Centrifuges



C - 852



R - 303



C - 854 / 8

Medico / Doctor Centrifuges

These are low priced units meeting the majority of the small capacity centrifuging requirements of medical practitioners, pathological laboratories and hospitals for routine centrifuging tests.

The units are fitted with resiliently mounted motor for vibration free performance having built in 5 step speed regulator. Available with Swing out or Angle rotor.

Technical Specifications

Parameters	C-852	C-854/4	C-854/6	C-854/8	R-303	R-304
Capacity (ml)	4 x15	4 x15	6 x15	8 x15	8 x15	
Type of Rotor	Swing Out				Angle	
No. of Tubes	4	4	6	8	8	
Max. Speed (RPM)	3500	3500	3500	3500	3800	
Max. RCF ('g')	1750	1850	1600	1600	1450	
W x D x H (mm)	295 x 295 x 250	310 x 310 x 295			210 x 210 x 250	
Supply: 220-240 Volts, 50 Hz Single Phase.						

Discontinued