



Neya 5C



Neya 5



NEYA 5/5C ROTOR AND ACCESSORIES


 HYBRID ROTOR	Aluminum Carriers (Set of 4) suitable for Hybrid Rotor			
				
MODEL	FA 50	SO 50	FA 15	SO 15
TYPE	FIXED ANGLE	SWING OUT	FIXED ANGLE	SWING OUT
CAPACITY PER CARRIER (ML)	1 X 50	1 X 50	2 X 15	2 X 15
TOTAL CAPACITY	4 X 50ML		8 X 15ML	
MAX. SPEED (RPM)	4700	4700	4700	4700
MAX. RCF (g)	3350	3560	3350	3560

IMAGE OF TUBE								
TUBE CAPACITY (ML)	50	50	15	15	15	10	5/7	5/7
TYPE OF TUBE	CONICAL	GLASS/PLASTIC	CONICAL	CONICAL	GLASS/PLASTIC	BLOOD COLLECTION TUBE	BLOOD COLLECTION TUBE (LONG)	BLOOD COLLECTION TUBE (SHORT)
CARRIER NAME	FA 50 SO 50	FA 50 SO 50	FA 50 SO 50	FA 15 SO 15	FA 15 SO 15	FA 15 SO 15	FA 15 SO 15	FA 15 SO 15
INSERT OR CUSHION NO.	C 50F	C 50R INCLUDED	RE 50-15F	C 15F	C 15R INCLUDED	C 10R	RE 15-5/7LB	RE 15-5/7SB
NO OF TUBES/CARRIER	1	1	1	2	2	2	2	2
NO OF TUBES/ROTOR	4	4	4	8	8	8	8	8
Ø X L (MM)	29 X 116	30 X 110 30 X 100	16.5 X 120	16.5 X 120	17 X 110 17 X 102	15.5 X 106	12.5 X 106	12.5 X 81
TUBE SUPPLIED	NO	YES	NO	NO	YES	NO	NO	NO