



RFV 340 Plus

These Deep Freezers are designed for quick freezing and storing of Blood Components, Serum, Vaccines, Biological & Medical Specimens, Clinical Samples, etc. at low temperatures.

Engineered for just a little bit more, the new REMI series of Upright Freezers gives you up to -25°C instead of the conventional -20° C.

Salient Features :

- Plastic inner chamber; Outer body pre painted G.I.
- Quick freezing function
- Defogger for removing condensation from Glass door.
- Solid state temperature controller
- Adjustable shelves
- Low energy consumption and noise levels
- Auto door closure
- Door open alarm
- Wheels provided with leveling jack.
- Integrated door handle
- Forced air circulation for Better uniformity within chamber
- In Built Voltage Stabilizer

Technical Data		
Model	RFV 245 Plus	RFV 340 Plus
Capacity (liters)	245	340
Temperature Setting Range	Ambient to -25 °C	
Door	1	1
No. of Trays	5	5
External dimensions W x D x H (mm)	625 x 690 x 1830	625 x 690 x 2085
Supply	220-240 V, 50 Hz, Single Phase	