






Rotor Heads & Accessories

Swing out Rotor (except SN 4-100 & SN 2-96 MP) & Angle Rotors are supplied with suitable SS tube as carriers

Swing Out Rotors for NEYA 6 & NEYA 4


Image of Tube	Metal Buckets suitable for Swing Out Rotor SN 4-100								
	B 100			B 9-5			B 9-5L		
									
Image of Tube									
Tube Capacity (ml)	100	50	50	15	15	10	5/7	5/7	5/7
Type of Tube	Glass	Glass	Conical	Glass	Conical	Blood Collection Tube	Blood Collection Tube (Long)	Blood Collection Tube (Short)	Blood Collection Tube (Long)
Bucket No.	B 100	B 100	B 100	B 100	B 100	B 100	B 100	B 9-5	B 9-5L
Insert or Cushion No.	C 100R Included	RE 100-50 (OLD - RA-157)	RE 100-50FR (OLD - RA-161F)	RE 100-15R(3) (OLD - RA-166)	RE 100-15FR (OLD - RA-160F)	RE 100-10(3) (OLD - RA-168)	RE 100-5/7L(4) (OLD - RA-167)	-	-
No of tubes/Insert	1	1	1	3	1	3	4	9	9
No of tubes/rotor	4	4	4	12	4	12	16	36	36
Ø x L (mm)	45.5 x 104	30 x 110	29 x 116	17 x 110	16.5 x 120	15.5 x 106	12.5 x 106	12.5 x 81	12.5 x 106
Tubes supplied	Yes	Yes	No	Yes	No	No	No	No	No

Swing out Rotor SN 4-100 with max. speed of 4500rpm, Max rcf (g) 3170 also available on Request

Swing Out Rotor SN 8-15 Max. Speed RPM 4500 Max RCF (g) 3440	Image of Tube							
	Tube Capacity (ml)	15	15	8	5/7	5/7	5	3
	Type of Tube	Glass	Conical	Glass	Blood Collection Tube (Long)	Blood Collection Tube (Short)	Glass	Glass
	Insert or Cushion No.	C 15R Included	C 15F	RE 15-8R (OLD - RA-154)	RE 15-5/7LR (OLD - RA-163)	RE 15-5/7SR (OLD - RA-162)	RE 15-5R (OLD - RA-153)	RE 15-3R (OLD - RA-152)
	No of tubes/Insert	1	1	1	1	1	1	1
	No of tubes/rotor	8	8	8	8	8	8	8
	Ø x L (mm)	17 x 110	16.5 x 120	12 x 100	12.5 x 106	12.5 x 81	12 x 75	10 x 75
	Tubes supplied	Yes	No	Yes	No	No	Yes	Yes


Swing Out Rotor SN 6-50 Max. Speed RPM 4500 Max RCF (g) 3530	Image of Tube					
	Tube Capacity (ml)	50	50	15	15	10
	Type of Tube	Glass	Conical	Conical	Glass	Blood Collection Tube
	Insert or Cushion No.	C 50R Included	C 50F	RE 50-15FR (OLD - RA-501F)	RE 50-15RR (OLD - RA-501)	RE 50-10R (OLD - RA-164)
	No of tubes/Insert	1	1	1	1	1
	No of tubes/rotor	6	6	6	6	6
	Ø x L (mm)	30 x 110	29 x 116	16.5 x 120	17 x 110	15.5 x 106
	Tubes supplied	Yes	No	No	Yes	No


Swing Out Rotor SN 24-15 Max. Speed RPM 4500 Max RCF (g) 3330	Image of Tube						
	Tube Capacity (ml)	15	8	5/7	5/7	5	3
	Type of Tube	Glass	Glass	Blood Collection Tube (Long)	Blood Collection Tube (Short)	Glass	Glass
	Insert or Cushion No.	C 15R Included	RE 15-8R (OLD - RA-154)	RE 15-5/7LR (OLD - RA-163)	RE 15-5/7SR (OLD - RA-162)	RE 15-5R (OLD - RA-153)	RE 15-3R (OLD - RA-152)
	No of tubes/Insert	1	1	1	1	1	1
	No of tubes/rotor	24	24	24	24	24	24
	Ø x L (mm)	17 x 110	12 x 100	12.5 x 106	12.5 x 81	12 x 75	10 x 75
	Tubes supplied	Yes	Yes	No	No	Yes	Yes


Swing Out Rotor SN 16-15 Max. Speed RPM 4500 Max RCF (g) 3400	Image of Tube							
	Tube Capacity (ml)	15	15	8	5/7	5/7	5	3
	Type of Tube	Glass	Conical	Glass	Blood Collection Tube (Long)	Blood Collection Tube (Short)	Glass	Glass
	Insert or Cushion No.	C 15R Included	C 15F	RE 15-8R (OLD - RA-154)	RE 15-5/7LR (OLD - RA-163)	RE 15-5/7SR (OLD - RA-162)	RE 15-5R (OLD - RA-153)	RE 15-3R (OLD - RA-152)
	No of tubes/Insert	1	1	1	1	1	1	1
	No of tubes/rotor	16	16	16	16	16	16	16
	Ø x L (mm)	17 x 110	16.5 x 120	12 x 100	12.5 x 106	12.5 x 81	12 x 75	10 x 75
	Tubes supplied	Yes	No	Yes	No	No	Yes	Yes




Rotor Heads & Accessories


Swing Out Rotor S 24-5 Max. Speed RPM 4400 Max RCF (g) 3200 	Image of Tube		
	Tube Capacity (ml)	5/7	5/7
	Type of Tube	Blood Collection Tube (Long)	Blood Collection Tube (Short)
	Insert or Cushion No.	Included	Included
	No of tubes/Insert	1	1
	No of tubes/rotor	24	24
	Ø x L (mm)	12.5 x 106	12.5 x 81
	Tubes supplied	No	No


Swing Out Rotor S 24-10 Max. Speed RPM 4400 Max RCF (g) 3200 	Image of Tube	
	Tube Capacity (ml)	10
	Type of Tube	Blood Collection Tube
	Insert or Cushion No.	No
	No of tubes/Insert	-
	No of tubes/rotor	24
	Ø x L (mm)	15.5 x 106
	Tubes supplied	No

Swing Out Rotor SN 2-96 MP Max. Speed RPM 4500 Max RCF (g) 2130 	Image of Tube	
	Tube Capacity (ml)	2 X 96 wells
	Type of Tube	Micro
	No of tubes/rotor	2
	W x D x L (mm)	127 x 84 x 18
	Tubes supplied	No

Angle Rotors for NEYA 6 & NEYA 4

Angle Rotor AN 4-100 Max. Speed RPM 6000 Max RCF (g) 4750 	Image of Tube							
	Tube Capacity (ml)	100	50	50	15	15	10	5/7
	Type of Tube	Plastic	Conical	Glass	Glass	Conical	Blood Collection Tube	Blood Collection Tube (Long)
	Insert or Cushion No.	C 100R Included	RE 100-50FR (OLD - RA-161F)	RE 100-50RR (OLD - RA-161)	RE 100-15R(3) (OLD - RA-166)	RE 100-15FR (OLD - RA-160F)	RE 100-10(3) (OLD - RA-168)	RE 100-5/7L(4) (OLD - RA-167)
	No of tubes/Insert	1	1	1	3	1	3	4
	No of tubes/rotor	4	4	4	12	4	12	16
	Ø x L (mm)	45.5 x 104	29 x 116	30 x 110	17 x 110	16.5 x 120	15.5 x 106	12.5 x 106
	Tubes supplied	Yes	No	Yes	Yes	No	No	No

Angle Rotor AN 6-50 Max. Speed RPM 6000 Max RCF (g) 4870 	Image of Tube					
	Tube Capacity (ml)	50	50	15	15	10
	Type of Tube	Plastic	Conical	Glass	Conical	Blood Collection Tube
	Insert or Cushion No.	C 50R Included	C 50F	RE 50-15RR (OLD - RA-501)	RE 50-15FR (OLD - RA-501F)	RE 50-10R (OLD - RA-164)
	No of tubes/Insert	1	1	1	1	1
	No of tubes/rotor	6	6	6	6	6
	Ø x L (mm)	30 x 100	29 x 116	17 x 110	16.5 x 120	15.5 x 106
	Tubes supplied	Yes	No	Yes	No	No

Angle Rotor AN 24-15 Max. Speed RPM 6000 Max RCF (g) 5070 	Image of Tube							
	Tube Capacity (ml)	15	15	8	5/7	5/7	5	3
	Type of Tube	Plastic	Conical	Glass	Blood Collection Tube (Long)	Blood Collection Tube (Short)	Glass	Glass
	Insert or Cushion No.	C 15R Included	C 15F	RE 15-8R (OLD - RA-154)	RE 15-5/7LR (OLD - RA-163)	RE 15-5/7SR (OLD - RA-162)	RE 15-5R (OLD - RA-153)	RE 15-3R (OLD - RA-152)
	No of tubes/Insert	1	1	1	1	1	1	1
	No of tubes/rotor	24	24	24	24	24	24	24
	Ø x L (mm)	17 x 102	16.5 x 120	12 x 100	12.5 x 106	12.5 x 81	12 x 75	10 x 75
	Tubes supplied	Yes	No	Yes	No	No	Yes	Yes