



Benchtop Tube Sealer WeldOn-B



WeldOn-B

WeldOn-B is a robust benchtop sealer with consistent performance. Ideally suited for blood component separation labs

WeldOn-B helps to keep the blood bag system completely sterile. Ensures quality sealing with clear rupture line for easy tear-apart segments.

Safety features & monitoring system makes it an ideal choice of users.



**Seal with clear break line
for snap-apart separation**



Splash Gaurd

Salient Features

- Senses the blood bag tube for automatic sealing process
- Anti-spark and overheat protection system to avoid damage to system & tubing
- Continuous status updates on LCD :
 - Ready to seal,
 - Sealing in process
 - Sealing Completed
- Easy to remove sealing head cap for quick cleaning
- Seal with clear break line for snap-apart separation
- Splash guard for technician safety
- LED Indication for cover opening

Technical Specifications

Parameters	WeldOn-B
Dimension (W x D x H) (mm)	220 x 200 x 120
Weight (Kg)	4
RF Frequency	40.68 MHZ
RF Output	100 watt
Seal Time	1.5 Sec (auto adjustable)
Sealing Type	Automatic

Input Voltage: 230V, 50 Hz